

No.1 Breathing Exercise & Breathing Monitoring

A global digital healthcare company that helps you breathe comfortably



Website



Amazon



Youtube

[주]지에이치이노텍

BiG
breathe

GHIT
Good Health Innotech

Tel | +82-51-711-4703 E-Mail | ghit@ghit.co.kr Web | www.ghit.co.kr

Address | Main Office 2F, B2, 68 Busandaehak-ro 50beon-gil, Geumjeong-gu, Busan (Jangjeon-dong, AVEC)

Branch Office 3F, 356 Bansong-ro, Geumjeong-gu, Busan

Taking the first step towards comfortable breathing, GH Innotek is here with you.

GH Innotek specializes in manufacturing and distributing respiratory products, reducing reliance on imported goods in the medical and sports sectors.

Our products are widely used in hospitals, welfare centers, and public health institutions, where they play a vital role in respiratory enhancement and rehabilitation programs aimed at improving public health.

Looking ahead, we are committed to developing a comprehensive respiratory care platform that will support remote medical services.

Under the brand name "Big Breathe," we are expanding our presence in the U.S., Europe, and Japan, achieving increased sales and solidifying our market position globally.

GH Innotek is dedicated to continuously developing and delivering innovative products that enhance the quality of life for people around the world.



- 2024**
 - 07 Designation of Outstanding R&D Innovation Products (O-PEP/IMT), Public Procurement Service
 - 06 Designated as a promising company (G-PASS company) entering the overseas procurement market, Public Procurement Service
 - ESG Leading Company Designation, Busan City
 - 01 Certification of Leading Strategic Industries (BioHealth Industry) , Busan City
- 2023**
 - 12 Award for contribution to the development of the medical industry : Busan Metropolitan City Mayor's Award
 - MIP/MEP FCC, CE Electromagnetic Compatibility Certification
 - 11 \$0.36 million in investment Korea SMEs and Startups Agency
 - 07 Procurement Innovation Product Pilot Purchase Supply (\$0.25 million)
 - 01 Selected as a procurement company in Massachusetts, U.S
- 2021**
 - 08 Venture Business Certification
 - 01 FDA Certification (5 items)
- 2020**
 - 05 ISO13485 Certification
- 2018**
 - 11 Corporate Research Institute Certification(KOITA)
 - 01 Research Institute Company Registration
- 2017**
 - 10 Establishment of GH Innotek Co.,Ltd



FDA Certification



ISO 13485:2016



FMDA Certification




US Trademark Registration Certificate




CE MDR DOC




Chinese patent Registration certificate




US Patent Registration Certificate



Venture Business Confirmation Certificate



Certificate of Innovative Product Designation



Medical Device Manufacturing License



Respiratory status monitoring device



Respiratory rehabilitation device



How to check FDA certification

- Scan the QR code on the left
- Enter GH Innotek in the Owner/Operator Name field



How to check EUDAMED registration

- Scan the QR code on the left
- Enter GH Innotek in the Manufacturer/Producer (and Authorized Representative) name field.

Medical Devices

Double action IMT/PEP



Product Specifications
19*4cm/50g
Key Features
3 Modes of Exercise

Product Description

- For COPD, respiratory patients, the elderly, and rehabilitation
- Strengthens respiratory muscles before and after surgery
- Prevents lung complications and clears airways
- Supports high-pressure training (up to 100cmH₂O) for sports

IMT Pressure Range	PEP Pressure Range
MAX 40cmH ₂ O	MAX 20cmH ₂ O

Double action V-PEP/IMT



Product Specifications
14*6cm/81g
Key Features
3 Modes of Exercise

Product Description

- Device for preventing and treating respiratory diseases like pneumonia
- Helps expel phlegm from the lungs/airways during exhalation
- Prevents chronic conditions from prolonged phlegm/coughing
- Four vibration modes, suitable for all ages

IMT Pressure Range / Vibrations per Second
7~46cmH ₂ O / Over 40Hz

Healthcare

Breath Care MIP/MEP



Product Specifications
4*5*9cm/46g
Key Features
Early detection of lung disease Web/app data monitoring

Product Description

- Self-diagnose respiratory function with muscle and lung capacity tests
- App integration for muscle age and lung capacity measurements
- Track exercises to improve lung function

MIP/MEP
Lung function test and breathing exercises in one!

Big Breathe S-IMT



Product Specifications
7.1*4.6cm/27g
Key Features
Inspiratory Muscle Strengthening

Product Description

- Adjustable pressure breathing trainer
- World's smallest resistance breathing device
- Slimming effect by strengthening the diaphragm
- Boosts endurance for better sports performance

1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
16	24	32	45	58	71	80	82	100

(Unit : cmH₂O)

Sports

Big Breathe IMT



Product Specifications
14*8cm/44g
Key Features
3 Types/Inspiratory Strength

Product Description

- 3-level pressure range for beginners to professional athletes
- Noticeable improvement in stamina and lung capacity
- Increased calorie burn based on pressure levels
- Precise pressure control for systematic training

모델	1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv	10Lv
Low	12	20	28	36	44	52	60	69	81	93
Middle	19	32	48	64	81	98	115	134	154	174
High	31	54	77	100	123	146	169	193	222	251

(Unit: cmH₂O)

Training Mask Brest Mask 1.0




Product Specifications
11*13*8.5cm/121g
Key Features
Boosts training performance Cool-down effect Maximizes calorie burn

Product Description

- Hands-free, usable with any activity
- Wide pressure range enhances calorie burn (Running Crew)
- Boosts lung capacity for improved endurance
- 28% more effective than airflow-based devices

1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
6	9	16	26	36	50	54	67	72

(Unit : cmH₂O)



Study Subjects

✓ 80 elderly participants (aged 71-74)
(Divided into 4 groups of 20)

- 1. Center-Based Training
- 2. Home-Based Training
- Group A : Big Breathe IMT/PEP
- Group C : Big Breathe IMT/PEP
- Group B : P Company IMT
- Group D : P Company IMT

Exercise Protocol

		Center-based Training		Home-based Training	
		Group A	Group B	Group C	Group D
Frequency	Supervised Exercise	2times/wk		2times/protocol	
	Self Exercise at home	2times/day For 8wks			
Intensity		40% of the initial MIP PEP : 20cmH ₂ O	40% of the initial MIP	40% of the initial MIP PEP : 20cmH ₂ O	40% of the initial MIP
Time		10 times/set, 20 sets (about 30min)			
Type		Combined IMTPEP	Philips IMT	Combined IMTPEP	Philips IMT

- **Center-Based Training**
 - Hospital visits twice a week for 8 weeks
- **Home-Based Training**
 - Supervised exercise sessions conducted during the first and last visits
- All groups perform home training for 8 weeks
 - One set of 10 repetitions / 10 sets, twice a day





Results

Outcomes	Group A n=17		Group B n=16		Group C n=19		Group D n=19	
	Change from baseline	P value	Change from baseline	P value	Change from baseline	P value	Change from baseline	P value
MIP cmH ₂ O	18.07 ± 17.87 (9.15 to 26.9)	0.001*	12.66 ± 12.76 (5.86 to 19.46)	0.001*	15.65 ± 15.96 (7.95 to 23.34)	0.000*	10.17 ± 9.85 (5.42 to 14.92)	0.000*
MIP % predicted	20.34 ± 20.35 (9.87 ± 30.80)	0.001*	14.04 ± 14.06 (6.55 ± 21.53)	0.001*	17.42 ± 17.82 (8.83 ± 26.01)	0.000*	11.20 ± 10.75 (4.54 ± 6.01)	0.000*
MEP cmH ₂ O	20.34 ± 20.35 (15.89 to 36.02)	0.000*	0.87 ± 21.66 (-10.66 to 12.41)	0.874	14.75 ± 21.58 (4.35 to 25.15)	0.008*	8.08 ± 22.57 (-2.79 to 18.97)	0.136
MEP % predicted	20.34 ± 20.35 (13.63 to 32.28)	0.000*	0.95 ± 19.57 (-9.47 to 11.38)	0.848	13.18 ± 19.00 (-4.02 to 22.34)	0.007*	6.96 ± 20.07 (-2.70 to 16.64)	0.148

Comparison of Before/After Effects of Using Big Breathe Product
 → MIP : **9.87~30.8% improvement**
 → MEP : **13.63~32.38% improvement**

Comparison of unidirectional breathing exercises (IMT, Inspiratory) with bidirectional exercises (Big Breathe)
 → MIP : **1.04~4.79% improvement**
 → MEP : **9.94~17.65% improvement**

MIP (Maximum Inspiratory Pressure) / MEP (Maximum Expiratory Pressure)



IMT Efficacy Summary - 1

Subjects

In a study of 200 post-COVID-19 patients, the IMT group outperformed the control group (standard breathing exercises).

- Dyspnea Severity: 28% Improvement
- FVC: 16% Improvement
- MIP: 16% Improvement

IMT Efficacy Summary - 2

Subjects

Among 59 post-COVID-19 patients, the IMT group showed better results than the control group (standard breathing exercises).

- Chest Expansion: 2x Improvement
- FEV1: 3x Improvement
- PEF: 3x Improvement

PEP Efficacy Summary - 1

Subjects

125 patients hospitalized due to COVID-19

- Experimental Method Improvement after use compared to before
- SpO2 / FiO2: 42.70% Increase

PEP Efficacy Summary - 2

Subjects

378 adults with COVID-19 symptoms

- Experimental Method PEP users outperformed the control group (standard breathing exercises)
- Results: Less chest tightness, less difficulty breathing, more energy, and higher daily activity levels

OPEP(VPEP) Efficacy Summary

Subjects

17 patients with non-cystic fibrosis bronchiectasis

- Experimental Method Pre/Post Comparison
- Results: Using OPEP devices helps prevent acute exacerbations in patients with bronchiectasis

	MIP-PRE cmH ₂ O	MIP-POST cmH ₂ O	MIP Increase %
Male	61.24	68.25	12.1
Female	41.25	50.38	22.12

	MEP-PRE cmH ₂ O	MEP-POST cmH ₂ O	MEP Increase %
Male	68.12	83.47	22.54
Female	53.13	67	26.12

Average increase compared to pre-use of IMT/PEP
MIP 7.96(14.5%) · MEP 14.88(23.5%) Increase



Development of a Simple Diagnostic Device for Respiratory Sarcopenia and Lung Capacity Using Pressure Sensors

Breath Care MIP/MEP

- ✓ Lung disease
 - ➔ 3rd leading cause of death, projected to be 1st by 2030
- ✓ Lung function is often found to have decreased by 30-40%
 - ➔ Early detection is crucial



BreathCare MIPMEP

Simple Diagnosis of Respiratory Muscle Weakness and Lung Capacity

Domestic Development of Testing Equipment (Price/Performance Advantage)

App Program

Self-diagnosis of respiratory pressure and lung capacity

Diagnosis of respiratory muscle lower limit

Includes a solution for diagnosing respiratory muscle age

Breathing exercises (strength training, endurance training)



2023 Innovative Product Pilot Purchase, Respiratory Exercise, and Function Measurement



Operation of the Respiratory Rehabilitation Program at Pusan National University Hospital
 Conducted at Pusan National University Hospital with education led by professional medical staff



Pusan National University Hospital (PNUH) Respiratory Training Device Trial
 Participants: Respiratory experts from PNUH
 Devices: Exercise bands, IMT LOW


*Average of three measurements used




Haeundae-gu Respiratory Rehabilitation Exercise Program
 Operated for 4 months with a 4-week session each month
 April, June, September, and November



Coverage of Sasang-gu Respiratory Rehabilitation Exercises
 Future expansion of respiratory rehabilitation exercise program to additional centers



Significant improvement in respiratory function, with increases in FVC, FEV1, MIP, MEP, and PEF, indicating enhanced respiratory muscle strength and capacity



PEP
Strengthening Expiratory Muscles for Secretion Clearance
 Helps remove phlegm from the lungs and airways
 Prevents lung complications

IMT + PEP
Simultaneous Use of IMT (Inhalation) and PEP (Exhalation)
 Reduces breathing exercise time compared to unidirectional devices

IMT
Strengthens the diaphragm by promoting abdominal breathing during respiration
 Delays lactic acid buildup → Improves exercise performance, Increases lung capacity

- Strengthens Respiratory Muscles
- Increases Diaphragm Thickness
- Improves Core Muscle Strength
- Increases Stroke Volume
- Enhances Lung Capacity

IMT Pressure Range	PEP Pressure Range
MAX 40cmH ₂ O	MAX 20cmH ₂ O

* IMT Pressure Setting: Recommended to be over 30% of MIP
 * Tolerance : Within ±2cmH₂O or ±10%

- ✓ Device for COPD, respiratory disease, and rehabilitation patients
- ✓ Strengthens respiratory muscles pre- and post-surgery
- ✓ Prevents lung complications and aids airway clearance
- ✓ Simultaneous inspiratory and expiratory muscle exercises

Features

<p>Technological Innovation</p> <p>Simultaneous Inspiratory and Expiratory muscle training</p>	<p>Convenient Cleaning</p> <p>Easily disassembled for cleaning</p>	<p>Price Competitiveness</p> <p>Dual-function device with one purchase</p>
---	---	---

Product Comparison

	Big Breathe	Other	Other
Breathing Method	Threshold Resistive	Flow resistive	Threshold Resistive
Exercise Efficiency	High	Low	High
Exercise Method	Bidirectional	Bidirectional	Bidirectional



PEP

- Helps prevent lung complications in patients after abdominal surgery
- Heart Rate : 6.8% reduction
- Respiratory Rate : 9% reduction

IMT

- FEV1 (Forced Expiratory Volume in 1 Second) : 7.9% increase
- FVC (Forced Vital Capacity) : 9.8% increase
- Pimax (Maximum Inspiratory Pressure) : 23.6% increase
- MVV (Maximum Voluntary Ventilation) : 22% increase
- Pm Peak, Pimax (Maximum Inspiratory Flow and Peak Pressure) : 26.6% increase
- MIP (Maximum Inspiratory Pressure) : 8% increase
- MEP (Maximum Expiratory Pressure) : 32% increase
- 6-Minute Walk Test (Distance covered in 6 minutes) : 4.9% increase

Resistance Exercise Efficacy

Threshold Resistive

- MVV 22%, MIP 29%, MEP 32%, YOYO Test 14%
 - 14-32% improvement in respiratory muscle strength
- Sustained respiratory function improvement

Threshold Resistive Group	VS	General Exercise Group
32% increase		Slight increase



ERS

RESPIRATORY EQUIPMENT AND DEVICES EXHIBITION MAGAZINE




Patient Groups	IMT Effects	PEP Effects
COPD, Pulmonary Emphysema	0	0
Asthma	0	
Heart Failure	0	0
NMD	0	0
Spinal Cord Injury	0	0
Hypertensi	0	0
Sleep Apnea	0	
Dysphagia	0	0
Vocal Cord Dysfunction	0	0
Lung Complication Prevention	0	
Spinal Cord Injury	0	
Low Back Pain	0	
Stroke	0	0

- ✓ Treats respiratory patients, the elderly, and those with COPD or in rehabilitation
- ✓ Helps remove phlegm during exhalation
- ✓ Lowers blood pressure for high blood pressure patients
- ✓ Strengthens lungs and improves lung capacity
- ✓ Prevents Pneumonia and Lung Complications
- ✓ Adjustable 4 ball weights suitable for all ages

Features

Technology Innovation Simultaneous use of vibration and inspiratory muscle strengthening (IMT and V-PEP combined)	Hygienic use Hygienic, attachable mouthpiece	Price competitiveness Dual-functionality with one product purchase Price competitiveness
---	--	---

Product Comparison

			
	Big Breathe	Other	Other
Breathing Method	Threshold resistive exhalation vibration type	Exhalation vibration type	Threshold resistive
Exercise Efficiency	Strengthens respiratory muscles / Chest wall vibration	Chest wall vibration	Strengthening the respiratory muscles
Exercise Method	Bidirectional	One-way	One-way



IMT
 Strengthens respiratory muscles by promoting inhalation and exhalation
 Delayed lactic acid release → Improved exercise performance, Increased lung capacity and endurance

IMT + V-PEP
 Simultaneous use of IMT (inhalation) and V-PEP (vibration during exhalation)
 Reduces respiratory training time compared to single-function products


V-PEP
 Vibrates to loosen phlegm for easier removal
 Supports airway clearance and prevents lung complications, Includes 4 interchangeable ball weights for all ages



Pressure Adjustment Range (1-9 Levels, 14 Steps)	IMT Pressure Range				Vibrations per second				
	1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
	7	11	16	21	26	31	36	40	46




(Unit : cmH₂O)
 * IMT Pressure Setting : Recommended to be over 30% of MIP (Maximal Inspiratory Pressure)
 * Tolerance : within ±2cm H₂O or ± 10%

Phlegm Removal	Increased Cardio-pulmonary Endurance	Improving Lung Capacity	Enhancing Exercise Endurance	Recovery From COVID-19
----------------	--------------------------------------	-------------------------	------------------------------	------------------------



- ✓ Adjustable pressure respiratory training device
- ✓ Strengthens diaphragm, reducing chest size
- ✓ Enhances cardiopulmonary endurance for sports and vocal performance
- ✓ **World's smallest pressure-resistance breathing device**

Product Comparison

			
	S-IMT	Other	Other
Breathing Method	Threshold resistive	Flow resistive	Flow resistive
Exercise Efficiency	High	Low	Low
Weight	27g	40g	139g

- Sports Training
- Diet
- Enhance Lung Capacity
- Respiratory Muscle Training
- Vocal Training



(Unit : cmH₂O)

S-IMT Pressure Adjustment Range (1-9 Levels, 14 Steps)								
1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
16	24	32	45	58	71	80	82	100

[Coming Soon] S-IMT PLUS Pressure Adjustment Range (1-9 Levels, 14 Steps)								
1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
25	37	59	88	122	143	166	182	200

* Tolerance : within ±2cm H2O or ± 10%

- ✓ Breathing exercise for respiratory muscle training
- ✓ Enhanced endurance for basic fitness and lung capacity
- ✓ Different levels for users: High intensity, Medium, and Low
- ✓ Increased calorie consumption based on pressure level

Features

Cardio Endurance


Improves cardiovascular Endurance for activities requiring significant breathing

Vocal Capacity

Helps improve singing ability by maximizing oxygen intake during training

Product Comparison

	Big Breathe	Other
Breathing Method	Threshold resistive	Threshold resistive
Exercise Efficiency	High	High
Product Weight	44g	72g
Anti-Infiltration Cap	Yes	No



(Unit : cmH₂O)



- Sports Training
- Diet
- Enhance Lung Capacity
- Respiratory Muscle Training
- Vocal Training

모델	1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv	10Lv
Low	12	20	28	36	44	52	60	69	81	93
Middel	19	32	48	64	81	98	115	134	154	174
High	31	54	77	100	123	146	169	193	222	251

* Tolerance : within ±2cm H₂O or ± 10%



- ✓ Increases exercise efficiency compared to airflow control
- ✓ Wide range of pressure adjustments strengthens respiratory muscles
- ✓ Enhances lung capacity and improves sports endurance
- ✓ Effective cool-down after exercise
- ✓ Increased calorie consumption according to pressure levels

Product Comparison

	BreST MASK	Other
Breathing Method	Threshold resistive	Flow resistive
Exercise Efficiency	Effective	Less
Input Pressure	Max 72 cmH ₂ O	7 cmH ₂ O

BreST MASK Wearing/Non-Wearing Group Comparison

*Rev Bras Med Esporte – Vol. 24, No 3 – Mai/Jun, 2018

	Wearing Group		Variables	Non-Wearing Group	
	PRE	POST		PRE	POST
	100±14	122±26	MVV %prev	102±13	100±25
Respiratory Muscle Strength	160±4	207±11	MIP cmH ₂ O	157±20	155±10
	151±44	200±10	MEP cmH ₂ O	153±20	153±35
YOYO TEST	660±45	750±21	DP meters	645±26	642±21

Threshold resistive
Major respiratory muscle strength improvement(14~32%)
 → Increase in respiratory muscle strength (sustained, not temporary)

Flow resistive
 No significant effect on respiratory muscle strength

(Unit : cmH₂O)

Pressure Adjustment Range 1~9 levels (14 stages total)								
1Lv	2Lv	3Lv	4Lv	5Lv	6Lv	7Lv	8Lv	9Lv
6	9	16	26	36	50	54	67	72

Breathing Training

Cool Down

Enhance Lung Capacity

High Intensity Exercise Efficiency






*MVV(Maximum voluntary ventilation) / MEP(Maximum expiratory pressure)
 MIP(Maximal inspiratory pressure) / DP(Distance traveled in YoYo Test)

- ✓ Self-diagnosis and medical use for respiratory muscles (MIP/MEP)
- ✓ Measures various inspiratory and expiratory lung capacities
- ✓ Checks respiratory function improvement using exercise records

Features

Respiratory Muscle Measurement	Lung Capacity Measurement	Breathing Exercises
Compare respiratory strength by age Measure respiratory muscle age Provides health risk warnings for respiratory function	Measurement of expiratory lung capacity [FVC(L), FEV1(L), FVC/FEV1(%), PEF(L/s)]	Muscle endurance Strength training

Product Comparison

			
	Breathe Care MIP/MEP	Other	Other
Price	\$150	\$1,100	\$7,400
Connection	Wireless (App/Web)	Wired	Wired
Size	Compact	Medium	Large
Respiratory Muscle Test	O	O	O
Lung Capacity Test	O	X	O
Remote Monitoring	O	X	X
Functions	Respiratory Function Test + Breathing Exercise, Remote Monitoring	Respiratory Muscle Measurement	Respiratory Function Measurement (Respiratory Muscles, Lung Capacity)



MIP/MEP Pressure Range

MAX $\pm 300\text{cmH}_2\text{O}$

- * MIP(Maximal inspiratory pressure)
- * MEP(Maximal expiratory pressure)

폐활량 측정 범위

8L 이상

- * MIP(Maximal inspiratory pressure)
- * MEP(Maximal expiratory pressure)

Increases Lung capacity

Strengthens Respiratory muscles

Family health monitoring

Early detection of lung disease

Measures Respiratory muscles and lung capacity

Key Features

Self-diagnosis Measures respiratory muscles and lung capacity
Breathing Exercises Improves endurance and muscle strength

Program

View data through the mobile app or web platform

Measurement Range

Respiratory Muscles : Within $\pm 300 \text{ cmH}_2\text{O}$
 Lung Capacity : Over 8L

Bluetooth Connectivity

Connects to the app via Bluetooth-enabled devices.

The screenshot displays three main sections of the app:

- Respiratory muscle measurement:** Shows MIP/MEP graphs and a table with values for 1st, 2nd, and 3rd breaths. Includes a 'Result Analysis' section with 'YES' and 'NO' options and a 'Pred' value of 92.
- Lung capacity measurement:** Shows Flow/Volume and Volume/Time graphs. Includes a 'Metrics' table with values for FVC, FEV1, FEV1/FVC, and PEF.
- Breathing exercises:** Shows progress for 'Breathing Strength Exercise' and 'Breathing Endurance Exercise' with circular progress indicators for accuracy and duration. Includes a 'Rest' section and a 'BREATHE IN' instruction.

The screenshot shows two tables from the web platform:

- Respiratory Muscle Results List:** A table with columns: Turn, Date, Name, MIP, MEP, and Select. It lists two entries for user 김희은.
- Lung Capacity Results List:** A table with columns: Turn, Date, Name, FVC, FEV1, and Select. It lists 11 entries for user 김희은.



MIP / MEP Measurement

- Compares average breaths by age and gender
- My respiratory muscle age
- Suggests hospital visit if below LLN*
- Confirms reproducibility**

Lung Capacity Measurement

Lung Function Test Index Measurement

1. FVC = Forced Vital Capacity
2. FEV1 = Forced Expiratory Volume in 1 Second
3. FEV1 / FVC = Standard indicator of airway obstruction
4. PEF = Peak Expiratory Flow rate

Breathing Exercise

- Respiratory muscle strength training
- Respiratory muscle endurance training

*LLN : Lower limit of normal

**Reproducibility: Determines accuracy by comparing values within 10% of the previous measurement.

Major business partners